

HA and β -TCP: Choosing a Phase Balance

RUIKANG Li (Neo)

Hongkong Kannry Industrial Co., Ltd. | www.kannrybio.com

Abstract. This note explains how to compare HA and β -TCP blends during material development. It includes a worked Ca/P calculation and the checks needed to distinguish composition from processing effects. Select a blend using measured phase composition and application requirements; Ca/P alone cannot predict resorption time.

Keywords: hydroxyapatite; β -TCP; phase balance; Ca/P; material development

1. Calculate the blend before comparing it

For ideal hydroxyapatite, $\text{Ca}_{10}(\text{PO}_4)_6(\text{OH})_2$, Ca/P is $10/6 = 1.667$. For tricalcium phosphate, $\text{Ca}_3(\text{PO}_4)_2$, it is $3/2 = 1.500$. A blend ratio expressed by mass cannot be inserted directly into an arithmetic average of these atomic ratios.

Worked example: combine 60 g stoichiometric HA and 40 g TCP. Using formula masses of approximately 1004.6 and 310.18 g/mol, calcium is $10(60/1004.6) + 3(40/310.18) = 0.984$ mol; phosphorus is $6(60/1004.6) + 2(40/310.18) = 0.616$ mol. The calculated Ca/P is about 1.597. This assumes pure, anhydrous phases and is not a measured batch result.

Use HA-only and β -TCP-only controls alongside a blend. A phase label alone cannot establish biological removal: a cell-excluding implantation study found different conversion behavior from a simple universal resorption ranking [1]. Compare the processed material, including its pore structure.

2. From phase balance to a trial

2.1 Set the comparison basis. Prepare a small matched series using one manufacturing route. State whether 60/40 means the weighed feed or the measured final phase ratio. If sintering changes the phases, assign the resulting material its measured identity; do not keep the feed ratio as the final specification.

2.2 Read two independent measurements. Pair quantitative XRD with bulk elemental analysis. A Ca/P value near the calculation supports elemental consistency but cannot rule out compensating secondary phases. A mismatch should trigger checks of moisture basis, carbonate or other substitutions, digestion recovery and the XRD phase model.

2.3 Choose an endpoint linked to the product. For a suspension, compare settling and redispersibility. For a ceramic body, compare open porosity and wet mechanical behavior. For an immersion study, hold medium volume, temperature and specimen geometry constant. Record phase changes and solution chemistry together; mass loss alone is insufficient.

3. Development decision

Advance the composition that meets the intended processing and physical requirements with the smallest unexplained variation. Select the phase ratio after this comparison, rather than assigning an arbitrary HA/TCP ratio to a desired clinical lifetime.

Table 1. Development comparison and interpretation controls.

Candidate	Purpose in the comparison	Control before interpretation
HA endpoint	Establish the HA-only response	Same size band and thermal history
β -TCP endpoint	Establish the TCP-only response	Confirm beta rather than alpha phase
60/40 mass blend	Test the intermediate composition	Verify phase fractions after processing
Processed blend	Separate blend from process effects	Measure porosity and exposed area

Scope. The 60/40 calculation is an idealized mass balance. It is neither a recommended implant composition nor evidence of in vivo performance.

References

1. Sakai S, Anada T, Tsuchiya K, Yamazaki H, Margolis HC, Suzuki O. Comparative study on the resorbability and dissolution behavior of octacalcium phosphate, β -tricalcium phosphate, and hydroxyapatite under physiological conditions. Dent Mater J. 2016;35(2):216-224. doi:10.4012/dmj.2015-255.

Prepared by: RUIKANG Li (Neo)



Reviewed by: Alex | Technical content: Approved

Review date: